


PROSTHETIC MANUAL

for Izenimplant System



PROSTHETIC MANUAL

for Izenimplant System

CONTENTS

Restoration Introduction

- 6 ZENEX IMPLANT SYSTEM
- 9 ZENEX System Overview
- 10 Prosthetic Type
- 11 Impression Type
- 12 Tightening Torque
- 13 Prosthetic Guide
- 14 Component & Instrument
- 16 Right Connection Checking Guide



21p

Cemented
Abutment

- 23 **Step 1** Separation of Cover Screw or Healing Abutment
- 24 **Step 2** Connect the Impression Coping
- 26 **Step 3** Impression Taking & Connect the Lab Analog
- 28 **Step 4** Working Model Production
- 29 **Step 5** Wax-Up, Casting & Porcelain Build-Up
- 30 **Step 6** Transfer Jig Production
- 31 **Step 7** Fastening of intraoral Abutment & installation of prosthesis



33p

Angled
Abutment

- 35 **Step 1** Separation of Cover Screw or Healing Abutment
- 35 **Step 2** Connect the Impression Coping
- 36 **Step 3** Impression Taking & Connect the Lab Analog
- 37 **Step 4** Working Model Production
- 37 **Step 5** Wax-Up, Casting & Porcelain Build-Up
- 38 **Step 6** Transfer Jig Production
- 39 **Step 7** Fastening of Abutment in oral cavity & installation of prosthesis



41p

Multi Straight
& Multi Angled
Abutment

- 44 **Step 1** Separation of Cover Screw or Healing Abutment
- 45 **Step 2** Connect the Multi Straight & Multi Angled Abutment in the oral cavity
- 46 **Step 3** Connect the Impression Coping
- 47 **Step 4** Impression Taking(Abutment level Impression taking)
- 48 **Step 5** Working Model Production
- 49 **Step 6** Wax-Up
- 50 **Step 7** Casting
- 51 **Step 8** Porcelain build up
- 52 **Step 9** Oxide film removal
- 53 **Step 10** Ceramic Crown Production
- 54 **Step 11** Delivering & Screwing



57p

Ball
Abutment

- 59 **Step 1** Separation of Cover Screw or Healing Abutment
- 60 **Step 2** Connect the Ball Abutment in the oral cavity
- 61 **Step 3** Impression Taking
- 62 **Step 4** Working Model Production
- 62 **Step 5** Wax Denture Production
- 63 **Step 6** Resin denture Production
- 65 **Step 7** Delivering


ZENEX IMPLANT SYSTEM

Designed for various types of bone

Post shape (Rounded top) optimized for digital dentistry system

Platform Switching

Two conical contact point






I system

11° tapered double contact connection



T system

Stage




Cover Screw

Select appropriate Cover Screw height upon depth of implant placement.

Select specification fits for fixture connection.

Tighten with 1.2 Hex Driver by hand.

Recommended tightening torque: 5~8Ncm




Healing Abutment

Use Healing Abutment fits for the diameter of abutment.

Use specification fits for fixture connection.

Tighten with 1.2 Hex Driver by hand.





Recommended tightening torque: 5~8Ncm



Healing Abutment	H	2.0	3.0	4.0	5.0	6.0	7.0	9.0
Abutment	G/H		1.0	2.0 or 3.0	2.0 or 3.0		Larger than 5.0	Larger than 5.0
Impression Coping	Type	Narrow	Narrow	Narrow	M / R / W	M / R / W	M / R / W	M / R / W

ZENEX System Overview

Single / Bridge Case

		Cemented Abutment	Angled Abutment	Multi Straight Abutment	Multi Angled Abutment
					
		2-Piece		3-Piece	
Prosthetic Type	Screw	✓	✓	x	x
	Cement	✓	✓	✓	✓
	Combination	✓	✓	✓	✓
Impression Type	Abutment Level	x	x	✓	✓
	Fixture Level	✓	✓	x	x

Single / Bridge Case

2 piece

Screw or cement or combination type rosthesis is possible with fixture level

Cemented / Angled Abutment

Screw or cement or combination type prosthesis is possible with fixture level impression, can be customized depending on oral environment and prosthesis type




3 piece

Multi Straight / Multi Angled Abutment

Cement or combination type prosthesis is possible with abutment level mpression and effective in bridge case with unfavorable path

* Contents above are general guideline from the company and products must be selected in consideration of oral environment, habits, fixture placement condition, clinical experience and aftermath.

Overdenture Case

		Ball Abutment	Multi Straight Abutment	Multi Angled Abutment
				
Prosthetic Type	Retentive Anchor	✓	✓	✓
	Bar Frame	✓	✓	✓
Impression Type	Abutment Level	✓	✓	✓
	Fixture Level	x	x	x

Overdenture Case

3 piece

Multi Straight / Multi Angled Abutment

Effective in the fabrication of overdenture using bar frame in abutment level impression

* Contents above are general guideline from the company and products must be selected in consideration of oral environment, habits, fixture placement condition, clinical experience and aftermath.

Prosthetic Type

Screw


Combined with abutment through casting and firing in fabrication process



- Screw hole is exposed above occlusal surface, therefore esthetics and occlusion have to be considered
- Prosthesis can easily be removed with screw, therefore there is no side effects from cement
- Errors can occur in bridge fabrication in casting or firing process
- Setting is affected severely by the fixture angle and adjacent teeth

Cement


Casted or fired separately from abutment in the fabrication process, and combined by cement



- There is no screw hole, therefore esthetic surface can be created
- Difficult to remove prosthesis
- Cement is difficult to remove and has chances for inflammation
- Passive fit in bridge is easy
- Relatively easy setting, only affected by adjacent teeth

Combination

· Casted or fired separately from abutment in the fabrication process, and combined by cement (same as cement type)




- Screw hole is exposed above occlusal surface, therefore esthetics and occlusion have to be considered
- Maintenance is easy because prosthesis can easily be removed with screw
- After connecting prosthesis with cement, cement can be removed completely outside the mouth, so there is no side effect from cement
- Passive fit in bridge is easy
- Setting is affected by the fixture angle and adjacent teeth but relatively easy compared to screw type

Impression Type


Abutment Level Impression

- Similar impression taking as natural teeth
- Bring abutment shape/position to working model (Impression taking is based on abutment information)
- Prosthetic process is relatively easy and convenient
- Close tray (ready made / stock tray) used
- Exclusive impression coping for each abutment is recommended




Fixture Level Impression Pick-up Type

- Bring fixture's connection/position to working model (impression taking is based on fixture information)
- Impression taking is relatively complicated but accuracy is better than transfer type
- Impression coping moves as one body with impression body
- Open tray (custom / individual tray) used



Fixture Level Impression Transfer Type

- Bring fixture's connection/position to working model (impression taking is based on fixture information)
- Convenient in posterior area with limited mouth opening
- Impression coping moves separately from impression body
- Close tray (ready made / stock tray) used



Tightening Torque

Recommended to use the tightening torque below

(Need regular maintenance for the abrasion, damage and functionality of components such as driver, torque wrench etc)

5 ~ 8Ncm



Cover Screw



Healing Abutment



Pick-Up Impression Coping

D432520



Transfer Impression Coping

D432520



Multi Pick-Up Impression Coping

D432520



Multi Transfer Impression Coping

D432520



Multi Healing Cap

D432520

20Ncm



Temporary Abutment



Multi Temporary Cylinder



Multi CCM Cast Cylinder



Multi Plastic Cast Cylinder

D432520



Multi Ti Link Cylinder

30Ncm



Cemented Abutment



Angled Abutment



Multi Straight Abutment




Multi Angled Abutment



Ball Abutment

Prosthetic Guide

Abutment Specification Selection




Guide Tip.

Emergence Profile Formation Tip

- Pre surgery planning is important since fixture depth decides abutment's G/H and H
- It is important to select abutment diameter similar to natural tooth's cervical area


Abutment Diameter Selection

● Diameter in cervical area
● Abutment Diameter




※ Natural teeth cervical area mesial-distal / buccal-lingual:
Based on smaller specification among standard specification

- When appropriate abutment specification for restoration was not selected impossible to create natural prosthesis contour like beside



ZENEX System Fixture
D \varnothing 4.5 / L 11.5mm

Cemented Abutment
D \varnothing 6.5 / P/H 5.5mm / G/H 2.0mm




ZENEX System Fixture
D \varnothing 4.5 / L 11.5mm

Cemented Abutment
D \varnothing 5.7 / P/H 5.5mm / G/H 2.0mm


Listing Number D432520

Component & Instrument


Prosthetic KIT (PDK)



Torque Wrench



Driver




※ Normally, perform rough connection by hand first and tighten in final torque with torque wrench

Right Connection Checking Guide

Cover Screw

- Misconnection happens by the bone near fixture or adjacent tissue and foreign substance
- Check right connection after removing interfering area with bone profiler



Healing Abutment


- If healing abutment and fixture has right connection, there is sealing on the top of taper area inside
- Misconnection happens by the bone near fixture or adjacent tissue and foreign substance
- Fixture failure can happen with plague and bacteria proliferation in gap
- Check right connection after removing interfering area with bone profiler



Impression Coping

Pick-up Impression Coping

- Misconnection occurs when fixture and hex do not connect accurately, or with interference from bone and tissue around fixture
- Check right connection by checking if coping body notch(A) matches with top of fixture or if there is gap inside the 11° taper area



Transfer Impression Coping

- Check right connection by checking if coping body notch(A) matches with top of fixture or if there is gap inside the 11° taper area


※ Transfer impression coping :

Guide pin will not be connected without accurately setting the hex, therefore reduce errors from users




Abutment


- Misconnection occurs when fixture and hex do not connect accurately, or with interference from bone and tissue around fixture
- Modify wrong hex setting with x-ray or use Bone profiler to remove interfering area and check right connection



Prosthetic Flow Chart




CEMENTED ABUTMENT



Prosthetic Process

- 23 **Step 1** Separation of Cover Screw or Healing Abutment
- 24 **Step 2** Connect the Impression Coping
- 26 **Step 3** Impression Taking & Connect the Lab Analog
- 28 **Step 4** Working Model Production
- 29 **Step 5** Wax-Up, Casting & Porcelain Build-Up
- 30 **Step 6** Transfer Jig Production
- 31 **Step 7** Fastening of intraoral Abutment & installation of prosthesis

Cemented Abutment



Mini, Regular & Wide

Abutment Features

The top part of the post is rounded, making it easy to fasten the zirconia crown.



Abutment for manufacturing Cement/Combination-retained type prosthesis

Select specification fits for fixture connection.

Customized by grinding

Fixture Level Impression


Tighten with 1.2 Hex Driver

Recommended tightening torque


Mini, Regular, Wide: 30Ncm

Abutment Diameter Selection

- Ø 4.0
- Ø 4.5
- Ø 5.2
- Ø 5.7
- Ø 6.5



Fixture Level Impression



Step 1



Separation of Cover Screw or Healing Abutment

Remove Cover Screw or Healing Abutment with 1.2 Hex Driver. At this time, connect dental floss to the spinner in the handle part of the Driver so that the Driver does not pass to the patient's neck. And prepare Impression Coping for the connection of the Fixture.

Healing Abutment



Cover Screw



Chair Side Step 1

Removal of Cover Screw or Healing Abutment

NOTE: Connect dental floss to the spinner in the handle part of the Driver so that the Driver does not pass to the patient's neck.

Step 2



Connect the Impression Coping

Transfer Impression Coping

Using 1.2 Hex Driver, connect the Transfer Impression Coping that matches the Fixture with the inside of the Fixture, and connect the Guide Pin.

1.2 Hex Driver



Transfer Impression Coping



Check the exact contact between the Impression Coping and the Fixture with X-ray.



1.2 Hex Driver



Pick-up Impression Coping



Pick-up Impression Coping

Using 1.2 Hex Driver, connect Pick-up Impression Coping that matches the Fixture with the inside of the Fixture and connect the Guide Pin.



Check the exact contact between the Impression Coping and the Fixture with X-ray.



Chair Side Step 2

Connect the Impression Coping

- Select the appropriate type of Impression Coping that matched the Fixture with the inside of the Fixture
- Using 1.2 Hex Driver, connect the Impression Coping and connect the Guide Pin
- Check the exact contact between the Impression Coping and the Fixture with X-ray.

Step 3



Impression Taking & Connect the Lab Analog

Transfer Impression Coping






After injecting the impression material using an Injection syringe around the Coping, the tray filled with the impression material is placed in the oral cavity to obtain an accurate impression.

After removing blood marks and other residues deposited in the impression body, separate the intraoral Transfer Impression Coping and connect it to the Lab Analog.

After contacting the connected Coping and Analog with the Transfer Impression Coping in the impression body, confirm the cross section accurately and deliver it to the Lab.

Transfer Impression Coping

Lab Analog



Pick-up Impression Coping



Lab Analog







Pick-up Impression Coping

After forming a suitable hole so that the head of the Guide Pin can be exposed outside the prepared individual tray, try the tray first to see if the head of the Guide Pin is visible through the hole.

Rubber impression material is injected without gap around Pick-up Impression Coping, and impression is obtained by accurately positioning the tray coated with the impression material.

After the impression material is hardened, unfasten the Guide Pin to remove the tray from the oral cavity, check for abnormalities in the impression body, remove bloodstains and other residues, and then deliver it to the lab.



Lab Side


Step 4

Working Model
Production

After passing the body making impressions in the doctor's office after checking the correct fastening of the Lab Analog form an Artificial Gum around Analog, and injecting the anhydrite is produced by the Working Model.

Accurately fasten the Cemented Abutment that matches the Analog on the work model.


Cemented Abutment



Step 5

Wax-Up, Casting & Porcelain Build-Up

Use Pattern Resin to make Resin Cap, and after Wax-Up, make PFM prosthesis in the usual way.



Lab Side

Step 6 (optional)

Transfer Jig Production

Remove Artificial Gum from the working model and connect the Abutment accurately using 1.2 Hex Driver. Next, build the Pattern Resin to make the Transfer Jig.

When making a prosthesis by directly fastening the Abutment on the model, it is necessary to make a Transfer Jig to reproduce the position of the Abutment on the model as it is in the oral cavity.

In particular, in the case of using a non-hex type Abutment, there is no repositioning function, so it is necessary to make an accurate Transfer Jig using pattern resin when receiving regardless of single/bridge.



Step 7




Fastening of intraoral Abutment & installation of prosthesis

Using a Transfer Jig, place the Abutment into the oral cavity accurately and fasten the Abutment with a 1.2 Hex Driver.

The correct connection between the Abutment and the Fixture is confirmed by X-ray.

The final tightening is to 30Ncm (to be tightened according to the recommended tightening torque value guided by Abutment) using a 1.2 Hex Driver and a Torque Wrench.

After checking the passive fit of the prosthesis margin, proper contact with the adjacent teeth, and occlusion with the antagonist teeth, block-out the screw hole and cement the final prosthesis on the Abutment.



1.2 Hex Driver



Torque Wrench




Cemented Abutment



ANGLED ABUTMENT


Prosthetic Process

- 35 Step 1** Separation of Cover Screw or Healing Abutment
- 35 Step 2** Connect the Impression Coping
- 36 Step 3** Impression Taking & Connect the Lab Analog
- 37 Step 4** Working Model Production
- 37 Step 5** Wax-Up, Casting & Porcelain Build-Up
- 38 Step 6** Transfer Jig Production
- 39 Step 7** Fastening of Abutment in oral cavity & installation of prosthesis




Angled Abutment

Mini, Regular & Wide



Abutment Features

The top part of the post is rounded, making it easy to fasten the zirconia crown.



Abutment for manufacturing Cement/Combination-retained type prosthesis

Various types of Angle

- 15° & 25° for Mini, Regular and Wide Fixture [Ø 3.5 ~ Ø 7.0]

Select specification fits for fixture connection.

Fixture Level Impression

Can be positioned in 12 directions by selecting A or B type

Tighten with 1.2 Hex Driver

Recommended tightening torque

Mini, Regular, Wide: 30Ncm

Step 1



Separation of Cover Screw or Healing Abutment

Remove Cover Screw or Healing Abutment with 1.2 Hex Driver. At this time, connect dental floss to the spinner in the handle part of the Driver so that the Driver does not pass to the patient's neck. And prepare Impression Coping for the connection of the Fixture.

Healing Abutment



Cover Screw



Step 2



Connect the Impression Coping

Transfer Impression Coping

Using 1.2 Hex Driver, connect the Transfer Impression Coping that matches the Fixture with the inside of the Fixture, and connect the Guide Pin.

1.2 Hex Driver



Transfer Impression Coping



Check the exact contact between the impression Coping and the Fixture with X-ray.



Step 3



Impression Taking & Connect the Lab Analog

Lab Analog






Transfer Impression Coping

After injecting the impression material using an Injection syringe around the Coping, the tray filled with the impression material is placed in the oral cavity to obtain an accurate impression.

After removing blood marks and other residues deposited in the impression body, separate the intraoral Transfer Impression Coping and connect it to the Lab Analog.

After contacting the connected Coping and Analog with the Transfer Impression Coping in the impression body, confirm the cross section accurately and deliver it to the Lab.



Step 4

Working Model Production

After passing the body making impressions in the doctor's office after checking the correct fastening of the Lab Analog form an Artificial Gum around Analog, and injecting the anhydrite is produced by the Working Model.

Accurately fasten the Angled Abutment that matches the Analog on the work model.

Angled Abutment



Step 5

Wax-Up, Casting & Porcelain Build-Up

Use Pattern Resin to make Resin Cap, and after Wax-Up, make PFM prosthesis in the usual way.



Lab Side


Step 6 (optional)

Transfer Jig Production

Remove Artificial Gum from the working model and connect the Abutment accurately using 1.2 Hex Driver. Next, build the Pattern Resin to make the Transfer Jig.

When making a prosthesis by directly fastening the Abutment on the model, it is necessary to make a Transfer Jig to reproduce the position of the Abutment on the model as it is in the oral cavity.

In particular, in the case of using a non-hex type Abutment, there is no repositioning function, so it is necessary to make an accurate Transfer Jig using pattern resin when receiving regardless of single/bridge.



Step 7




Fastening of Abutment in oral cavity & installation of prosthesis

Using a Transfer Jig, place the Abutment into the oral cavity accurately and fasten the Abutment with a 1.2 Hex Driver.

The correct connection between the Abutment and the Fixture is confirmed by X-ray.

The final tightening is tightened to 30Ncm (to be tightened according to the recommended tightening torque value guided by Abutment) using a 1.2 Hex Driver and a Torque Wrench.

After checking the passive fit of the prosthesis margin, proper contact with the adjacent teeth, and occlusion with the antagonist teeth, block-out the screw hole and cement the final prosthesis on the Abutment.



1.2 Hex Driver




Torque Wrench



Angled Abutment



Prosthetic Flow Chart




MULTI STRAIGHT & MULTI ANGLED ABUTMENT

Prosthetic Process

- 44 **Step 1** Separation of Cover Screw or Healing Abutment
- 45 **Step 2** Connect the Multi Straight & Multi Angled Abutment in the oral cavity
- 46 **Step 3** Connect the Impression Coping
- 47 **Step 4** Impression Taking(Abutment level Impression taking)
- 48 **Step 5** Working Model Production
- 49 **Step 6** Wax-Up
- 50 **Step 7** Casting
- 51 **Step 8** Porcelain build up
- 52 **Step 9** Oxide film removal
- 53 **Step10** Ceramic Crown Production
- 54 **Step 11** Delivering & Screwing



Multi Straight Abutment



Mini, Regular & Wide

Abutment for manufacturing screw-retained prosthesis in Multiple Case


Same platform as Multi Angled Abutment

Move into internal oral part by using exclusive Abutment Carrier
(Code: MSACR48)

Tighten with exclusive driver (Code: MSADSR20)

Recommended tightening torque: 30Ncm

Multi Angled Abutment



Mini, Regular & Wide

Abutment for manufacturing screw-retained prosthesis in Multiple Case

Abutment of various angles (17°, 30°) for various angled of implant insertion path

Same platform as Multi Straight Abutment

Connect by using exclusive Abutment Carrier (Code: MAACRMC)

Tighten with 1.2 Hex Driver

Recommended tightening torque: 30Ncm

Multi Angled Abutment Screw (MAASSR23 for Mini, Regular and Wide) included

Step 1 

Separation of Cover Screw or Healing Abutment

Remove Cover Screw or Healing Abutment with 1.2 Hex Driver. At this time, connect dental floss to the spinner in the handle part of the Driver so that the Driver does not pass to the patient's neck. And prepare Impression Coping for the connection of the Fixture.

Healing Abutment



Cover Screw



Multi Straight



Multi Angled



Step 2



Connect the Multi Straight & Multi Angled Abutment in the oral cavity

After connecting the Multi Straight Abutment to the Fixture with the Multi Straight Abutment Driver, check the connection between the Abutment and the Fixture with X-ray, and tighten it with 30Ncm using a Torque Wrench.

Multi Straight



After connecting the Multi Angled Abutment to the Fixture with a 1.2 Hex Driver, check the connection between the Abutment and the Fixture with X-ray, and tighten the screw with 30Ncm using a Torque Wrench.

Multi Angled



Step 3

Connect the Impression Coping

Multi Transfer Impression Coping

Using 1.2 Hex Driver, connect Multi Transfer Impression Coping to Abutment.

Multi Transfer Impression Coping




Multi Straight



Multi Angled



Check the exact contact between the impression Coping and the Abutment with X-ray.



Step 4

Impression Taking

(Abutment level Impression taking)

Multi Lab Analog




Multi Transfer Impression Coping





After injecting the impression material using an Injection syringe around the Coping, the tray filled with the impression material is placed in the oral cavity to obtain an accurate impression.

After removing blood marks and other residues deposited in the impression body, separate the intraoral Multi Transfer Impression Coping and connect it to the Multi Lab Analog.




After contacting the connected Coping and Analog with the Transfer Impression Coping in the impression body, confirm the cross section accurately and deliver it to the Lab.



Multi Straight



Multi Angled



Lab Side

Step 5

Working Model
Production

Tighten the Multi Healing Cap to protect the Abutment until the prosthesis is installed.

Checking whether the Coping is well located in Multi Lab Analog.

After injecting artificial gum around the Analog, when it is hardened, pour stone to make a working model.

Multi Healing Cap



Connect the
Impression Coping
and Multi Lab Analog

Multi Lab Analog
located on the
impression body

Artificial Gum
Formation

Plaster
Injection after
boxing

Complete the
work model

Multi Straight



Multi Angled




Step 6

Wax-Up




Place the Multi CCM Cast Cylinder above the Abutment and fasten the cylinder screw with 20Ncm using a 1.2 Hex Driver.

After adjusting the plastic sleeve to the appropriate height, perform wax-up for the metal structure of the prosthesis.



Multi CCM Cast Cylinder



Multi Straight



Multi Angled



Lab Side

Step 7

Casting

The sprue is mounted on the margin, when forming a sprue for casting,

The Abutment metal part and the adjacent connection part are compensated with wax as much as possible.

It is recommended to use Ni-Cr alloy for casting metal.

Prohibition of use of Co-Cr alloy (excessive oxide film formation and casting shrinkage)

CCM Cast Abutment has different casting characteristics compared to Gold UCLA Abutment, so an oxide film is generated on the metal part after casting.

Multi Straight



Multi Angled



Step 8

Porcelain build up

Porcelain building up and firing on the casting

Polishing and polishing working in general

Checking if there is anything wrong with the working model


Removal of oxide film generated during casting and porcelain firing

- ① Block out with utility wax, etc., except for the metal part where the oxide film is generated.
- ② Primary removal of oxide film by blasting with a glass bead (4~6 bar) :
Do not use rubber wheel / point (damage to the connection part)

Multi Straight



Multi Angled



Lab Side

Step 9

Oxide film removal

- ① Remove the blocked out part : Final removal of oxide film by high polishing with rouge applied to cotton
- ② After high polishing, Ultrasonic or steam cleaning.

Multi Straight




Multi Angled



Step 10

**Ceramic Crown
Production**


The planned ceramic prosthesis is fabricated in the usual way.



Multi Straight



Multi Angled




Step 11

Delivering & Screwing

After checking the margin passive fit of the final prosthesis and checking occlusion and esthetics, first fasten with a 1.2 Hex Driver and a torque wrench with a 1.2 Hex Driver to completely fasten the prosthesis with 20Ncm.

After covering the protective material over the screw head, the access hole finishes the occlusal surface with resin in the oral cavity.



1.2 Hex Driver



Torque Wrench




Multi Straight



Multi Angled



Prosthetic Flow Chart




BALL ABUTMENT



Prosthetic Process

- 59 **Step 1** Separation of Cover Screw or Healing Abutment
- 60 **Step 2** Connect the Ball Abutment in the oral cavity
- 61 **Step 3** Impression Taking
- 62 **Step 4** Working Model Production
- 62 **Step 5** Wax Denture Production
- 63 **Step 6** Resin denture Production
- 65 **Step 7** Delivering

Ball Abutment



Mini, Regular & Wide

Abutment for overdenture using O-ring attachment

Compensation of mounting angle up to 20°

Tighten with exclusive Ball Abutment Driver (Code: BAD24)

Recommended tightening torque: 30Ncm

Step 1



Separation of Cover Screw or Healing Abutment

Remove Cover Screw or Healing Abutment with 1.2 Hex Driver. At this time, connect dental floss to the spinner in the handle part of the Driver so that the Driver does not pass to the patient's neck. And prepare Impression Coping for the connection of the Fixture.

Healing Abutment



Cover Screw



1.2 Hex Driver



Step 2




Connect the Ball Abutment in the oral cavity

Determine the proper height of the Ball Abutment by measuring the depth of the gingival tissue on the Fixture.

The shoulder of the Abutment should be positioned above the tissue (about 1.5-2mm).

Connect the Ball Abutment to the Fixture with the Driver for Ball Abutment.

After confirming by X-ray, connect the Torque Wrench to the Ball Abutment Driver and tighten it to 30Ncm. (※ When re-fastening the Healing Abutment after taking an impression, fasten it only with finger pressure.)



Ball Abutment



Ball Abutment Driver



Torque Wrench



Step 3



Impression Taking

The rubber impression material is first injected in detail around the Ball Abutment, and then the impression material is filled in the prepared individual tray and placed in the oral cavity to obtain an impression.

After checking for abnormalities in the impression body, send it to the lab.


If there is a denture currently installed, it can be used as a Temporary denture by modifying the part where the Ball Abutment is fastened.



Preliminary procedure: Before installing the Ball Abutment, take an impression of the edentulous extension with alginate impression material and send it to the workshop to make a personal tray (※ 2mm more relief than the height of the Abutment).

Scratch is formed on the border so that the impression material can be attached well.

Ball Abutment Lab Analog



Lab Side



Step 4

Working Model
Production

When the impression body is delivered for the pore process, the lab Analog is pushed into the inner surface of the impression body until it is completely inserted into the Ball portion.

Make a working model by carefully pouring stone so that the Analog position does not move.

Base plate and wax occlusal rim for occlusal acquisition are made and sent to the clinic with the model.



Lab Side

Step 5

Wax Denture Production

The occlusal rim is placed in the oral cavity to obtain an intermaxillary occlusion and sent back to the studio.

In the workshop, denture teeth are arranged on the wax rim according to the occlusal record sent.

It is sent back to the doctor's office to check the occlusion of the arranged teeth and check the functionality and aesthetics of the denture.

(※ If corrections are made, set up a new occlusal record and retry until fit is achieved.)






Step 6

Resin denture Production

When the oral fit for the wax denture is completed, the final resin denture is fabricated.

Check the black Lab O-ring in the Ball Abutment Retainer or Ball Abutment Cap (※ Make the smaller of the Ball Abutment Retainer openings be the occlusal side) and place it on the Lab Analog.

Block-out the lower part of the retainer with a putty to prevent the acrylic resin from flowing under the Ball Abutment Retainer, and make it about 2mm higher than the retainer to give mobility on the denture base.



The dentures are buried together with the Ball Abutment Retainer Assembly in place, and flasking, curing, and finishing are performed as usual to complete the fabrication.



Step 7

Delivering



Replace the Black Lab O-ring inside the Ball Abutment Retainer with the orange O-ring for final.

Adjust occlusal and tissue contact areas as needed.


While attaching attachments, inform the patient about oral hygiene and cautions when attaching and detaching dentures.

Replace O-rings when fatigue accumulates and cannot function.

It is replaced approximately once a year.



O-ring



IZ-PMA-02 REV.00 (JUN.22)



Izenimplant Co.,Ltd.
www.izenimplant.com

HQ & Factory

1, 2Dong, 26-32, Suworam 4-gil, Seotan-myeon,
Pyeongtaek-si, Gyeonggi-do, 17703, Republic of
Korea

Sales Office

#1013, #1014 Central-Tower II, 158, Magokseo-ro,
Gangseo-gu, Seoul, 07788
Republic of Korea

Tel. +82-31-662-0657 (HQ)
+82-2-3664-3666 (Sales)
Fax. +82-31-662-0658